

IN 1954 German scientist Reinhardt Voll, possibly by accident, noticed, that a patient, which touched the medication with his hand – at the moment of measuring his biologically active dots potential – had the measurement results changed. Thus, the invention was made – the medications, even at distance, could change the energetic condition of a person. In some time also the bio-resonance theory appeared. Bio-resonance therapy is distinctive due to the fact, that for treating the human no outer affections are performed. Using the special equipment pathological frequencies can be separated (from all frequencies range) and inverted 180 degrees. Human organism receiving such waves come back to the normal condition. Human is a multi-connection hierarchy system of direction, and the signals in itself are divided into harmonic ones, called physiologically healthy oscillations, and disharmonic ones, called pathological. The later ones distort the direction processes in the organism and may result in disease. Such oscillation should be eliminated in course of treatment. There are, actually, noise oscillations, nevertheless even if they disturb the direction processes, they do not result in pathology.

Bio-resonance therapy (BRT) is a new direction in the medical and prophylactic medicine. The device of bio-resonance inspection and frequency wave characteristics of different organs and systems work with the oscillations of a human organism. The electro-magnetic and wave oscillations of the patient are used as therapeutic agents. The Sensitiv Imago devices neutralize the pathological oscillations of the patient and activate one's energetic and informational, biochemical and physiological balance. The bio-resonance therapy is most effective in cases when to achieve the success in treatment by traditional methods is impossible or requires much time, or has some other difficulties.

Sensitiv technology is based on the well-known in medicine biomagnetic resonance phenomenon. In the scope of this theory, the human's organism like any other is considered as an open oscillatory circuit, which emits faint electromagnetic waves with the unique frequencies range. As these signals are ones of the electromagnetic nature, they can be registered using the electrodes or inductors and can be transmitted via the electric wires to the special devices.

Implementation of different outer substances in the organism or some pathological processes development distort the internal harmony and results in the changes of the frequencies range, especially, in the pathologic oscillations appearance. Now let's imagine, that a microbe or a toxin excreted by the parasites is constantly circulating in the organism and also has its own oscillation range. If a signal with analogy frequency characteristics, the phenomenon of resonance appears, which sensitive devices may register. Connected to the patient electrodes and inductors, the required frequencies range is transmitted. Using a special unit in the device and modern processors, the presence of different pathologic focuses, imbalance of the neural-reflector conductivity of tissues in human organism can be revealed with accuracy of nearly 97% through the non-invasive method. In order to make this method possible a great peace of mutual work was made before by the physicists, doctors, programmers. First of all, the electronic analogue of the resonance frequency should be created, for example, of the erosive gastritis, or rheumatoid arthritis, then to create an electronic version, to insert it in the electronic selector of the software, to analyze hundreds of analogue spectrographic indicators of different age groups of patients – and only after this – to determine highly precise analogue of this or that pathological process or of the bacterial-viral agent. In this case, the specific phenomenon is registered – scientists call that phenomenon “memory of the wave field”. Viruses, bacteria, fungi, protozoa, various parasites – all this microscopic armada cannot be detected in few minutes. Only the electronic analogs are needed. These electronic analogs are written into the software database after precise multi-level statistic researches and comparative analysis of the spectrographic data. Sensitiv Imago equipment provides the modulation of the low, medium and high frequencies – thus providing the whole spectrum of possibilities to perform both diagnostic and treatment effect on the patient. The thing is that low-frequency oscillations, when electrodes are contacting the patient, provide the impulse flow from one electrode to another – via the tissues. High-frequency impulse in turn is like “flowing over the upper layer, skin surface”. Patient's body surface in this case is considered mostly not a conductor, but a wave transmitter (like the coaxial cable, which is transferring the antenna-received signal to the TV or other receiver device). Thus, Sensitiv Technology, as a treatment method, can be considered not electropuncture, but the electro-magnetic, frequency-wave diagnostic and therapy. When the electric impulse, conducted by the electrodes, is a source of the secondary appearing electro-magnetic field, and the affection comes onto the patient, all pathological processes and stirrers it is aimed on.